

Activity: pH and Lemonade

Materials for Team of 4:

- lemonade concentrate, thawed
- distilled water
- pH paper
- pH color chart
- graduated cylinder
- beaker or large jar
- baby food jars

Procedures:

- Using distilled water, mix the lemonade as directed on the package.
- Determine the pH of the lemonade. Record this value. Determine and record the pH of distilled water.
- Measure 50 ml of the lemonade; add 50 ml of distilled water. Predict the pH of this solution. Test your prediction.
- Measure 25 ml of the lemonade; add 75 ml of distilled water. Predict and test the pH of this solution.
- Measure 10 ml of lemonade; add 90 ml of distilled water. Test and record the pH of this solution. Label this solution A.
- Take 10 ml of solution A and a different amount of distilled water of your own choosing. Predict and test the pH of your solution. Make 3 more combinations of distilled water. Start each trial with 10 ml of solution A. Test and record the pH. Graph the results. Compare your results with at least 2 other teams.
- After completing your testing, you may drink your lemonade solutions.

Interpretations:

- How does adding water influence pH?
- Predictions do not have to be exactly correct to be good. Did your predictions get better with practice? Why?
- Why do you think some lemonade tastes better than others?

Enrichment:

Write a poem or advertisement about pH and lemonade.